

**Certificate of Analysis
GM501 Mineral Oil**

Batchnumber	OIL-2112
Production date	01 June 2021
Expiry date	31 December 2022
Total volume (mL)	490000

	Results	Specifications
Appearance	PASSED	PASSED
Sterility test <small>Sterilised by sterile filtration over 0,2µm filter</small>	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	20	< 30 cP at 30°C
Specific gravity	0,849	< 0,86 g/ml
Endotoxin	< 0,1	< 0,1 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay <small>(1-cell to expanded blastocyst after 120 hours)</small>	95,24	≥ 80 %

Release date	18 June 2021
Released by	Vanessa Bierlich

Signature