

Certificate of Analysis GM501 Mineral Oil

Batchnumber	OIL-2211
Production date	05 April 2022
Expiry date	31 October 2023
Total volume (mL)	400000

	Results	Specifications
Appearance	PASSED	PASSED
Sterility test <small>Sterilised by sterile filtration over 0,2µm filter</small>	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	22,7	< 30 cP at 30°C
Specific gravity	0,8488	0,83-0,86 g/ml
Endotoxin	< 0,1	< 0,1 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay 1 <small>(1-cell to expanded blastocyst after 120 hours)</small>	100	≥ 80 %
Mouse Embryo Assay 2 <small>(1-cell to expanded blastocyst after 120 hours)</small>	100	≥ 80 %

Release date	21 April 2022
--------------	---------------

Released by	Dr. Julia Heinzmann
-------------	---------------------

Signature

p.p. Julia Heinzmann

REV 08_00 / 2021-09-15