

**Certificate of Analysis
GM501 Mineral Oil**

Batchnumber OIL-2217
 Production date 02 June 2022
 Expiry date 31 December 2023
 Total volume (mL) 191800

	Results	Specifications
Sterility test Sterilised by sterile filtration over 0,2µm filter	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	29	< 60 cP at 30°C
Specific gravity	0,85	0,83-0,86 g/ml
Endotoxin	< 0,1	< 0,1 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay 1 (1-cell to expanded blastocyst after 120 hours)	100	≥ 80 %
Mouse Embryo Assay 2 (1-cell to expanded blastocyst after 120 hours)	91	≥ 80 %

Release date 22 June 2022

Released by Vanessa Bierlich

Signature