

## Certificate of Analysis GM501 Mineral Oil


Batchnumber	OIL-2302
Production date	17 January 2023
Expiry date	31 January 2024
Total volume (mL)	254000

	Results	Specifications
Appearance	PASSED	PASSED
Sterility test <small>Sterilised by sterile filtration over 0,2µm filter</small>	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	21,5	≤ 30 cP at 30°C
Specific gravity	0,84	0,83-0,86 g/ml
Endotoxin	< 0,1	< 0,1 EU/ml
Peroxide value (POV)	< 0,003	< 0,1 mEq/kg
Mouse Embryo Assay 1 <small>(1-cell to expanded blastocyst after 120 hours)</small>	100	≥ 80 %
Mouse Embryo Assay 2 <small>(1-cell to expanded blastocyst after 120 hours)</small>	100	≥ 80 %

Release date 14 February 2023

Released by Nicole Geiger

Signature



REV 09\_00 / 2022-07-29