


**Certificate of Analysis
GM501 Mineral Oil**

Batchnumber OIL-2418
 Production date 27 May 2024
 Expiry date 30 November 2025
 Total volume (ml) 180.000

Parameter	Results	Specifications
Solid paraffin	PASSED	PASSED
Acidity	NEUTRAL	NEUTRAL
Readily carbonisable substances	PASSED	PASSED
Density	0,842	0,83-0,86 g/ml
Peroxide value measurement (POV)	< 0,06	< 0,1 mEq/kg
Viscosity at 30°C	20,0	< 30 cP
Viscosity at 37°C	15	-
Endotoxin	< 0,1	< 0,1 EU/ml
Sterility test (according to the current Ph.Eur.2.6.1 / USP <71>)	No growth	NO GROWTH
Mouse Embryo Assay 1 (1-cell to expanded blastocyst after 120 hours)	100	≥ 80 %
Mouse Embryo Assay 2 (1-cell to expanded blastocyst after 120 hours)	100	≥ 80 %

Release date 12 June 2024
 Released by Nicole Geiger
 Signature 

REV 10_00 / 2023-09-15

