

**Certificate of Analysis  
GM501 Mineral Oil**


Batchnumber OIL-2423  
 Production date 26 June 2024  
 Expiry date 31 December 2025  
 Total volume (ml) 180.000

Parameter	Results	Specifications
Solid paraffin	PASSED	PASSED
Acidity	NEUTRAL	NEUTRAL
Readily carbonisable substances	PASSED	PASSED
Density	0,84	0,83-0,86 g/ml
Peroxide value measurement (POV)	< 0,06	< 0,1 mEq/kg
Viscosity at 30°C	19,5	< 30 cP
Viscosity at 37°C	15	-
Endotoxin	< 0,1	< 0,1 EU/ml
Sterility test (according to the current Ph.Eur.2.6.1 / USP <71>)	NO GROWTH	NO GROWTH
Mouse Embryo Assay 1 (1-cell to expanded blastocyst after 120 hours)	100	≥ 80 %
Mouse Embryo Assay 2 (1-cell to expanded blastocyst after 120 hours)	100	≥ 80 %

Release date 11 July 2024

Released by Nicole Geiger

Signature



REV 10\_00 / 2023-09-15