


**Certificate of Analysis  
GM501 Mineral Oil**

Batchnumber OIL-2430  
 Production date 09 September 2024  
 Expiry date 31 March 2026  
 Total volume (ml) 180.000

Product was sterilized by aseptic processing techniques

<b>Parameter</b>	<b>Results</b>	<b>Specifications</b>
Solid paraffin	PASSED	PASSED
Acidity	NEUTRAL	NEUTRAL
Readily carbonisable substances	PASSED	PASSED
Density	0,8408	0,83-0,86 g/ml
Peroxide value measurement (POV)	< 0,06	< 0,1 mEq/kg
Viscosity at 30°C	19,9	< 30 cP
Viscosity at 37°C	15,2	-
Endotoxin	< 0,1	< 0,1 EU/ml
Sterility test (according to the current Ph.Eur.2.6.1 / USP <71>)	NO GROWTH	NO GROWTH
Mouse Embryo Assay 1 (1-cell to expanded blastocyst after 120 hours)	100	≥ 80 %
Mouse Embryo Assay 2 (1-cell to expanded blastocyst after 120 hours)	95,24	≥ 80 %

Release date 30 September 2024  
 Released by Nicole Geiger  
 Signature



REV 10\_00 / 2023-09-15