

Certificate of Analysis GM501 Mineral Oil

Batchnumber	OIL-2615
Production date	01 June 2026
Expiry date	31 December 2027
Total volume (ml)	260.560


Product was sterilized by aseptic processing techniques

Parameter	Results	Specifications
Solid paraffin	PASSED	PASSED
Acidity	NEUTRAL	NEUTRAL
Readily carbonisable substances	PASSED	PASSED
Density	0,85	0,83-0,86 g/ml
Peroxide value measurement (POV)	< 0,06	< 0,1 mEq/kg
Viscosity at 30°C	21,8	< 30 cP
Viscosity at 37°C	16,2	-
Endotoxin	< 0,1	< 0,1 EU/ml
Sterility test (according to the current Ph.Eur.2.6.1 / USP <71>)	No growth	NO GROWTH
Mouse Embryo Assay 1 (1-cell to expanded blastocyst after 120 hours)	90	≥ 80 %
Mouse Embryo Assay 2 (1-cell to expanded blastocyst after 120 hours)	100	≥ 80 %

Release date	18 June 2026
--------------	--------------

Released by	Nicole Geiger
-------------	---------------

Signature



REV 10_00 / 2023-09-15