

Certificate of Analysis
GM501 Mineral Oil

Batchnumber OIL-2621
Production date 27 July 2026
Expiry date 31 January 2028
Total volume (ml) 180.240

Product was sterilized by aseptic processing techniques

Parameter	Results	Specifications
Solid paraffin	PASSED	PASSED
Acidity	NEUTRAL	NEUTRAL
Readily carbonisable substances	PASSED	PASSED
Density	0,85	0,83-0,86 g/ml
Peroxide value measurement (POV)	< 0,06	< 0,1 mEq/kg
Viscosity at 30°C	21,8	< 30 cP
Viscosity at 37°C	16,9	-
Endotoxin	< 0,1	< 0,1 EU/ml
Sterility test (according to the current Ph.Eur.2.6.1 / USP <71>)	NO GROWTH	NO GROWTH
Mouse Embryo Assay 1 (1-cell to expanded blastocyst after 120 hours)	100	≥ 80 %
Mouse Embryo Assay 2 (1-cell to expanded blastocyst after 120 hours)	95	≥ 80 %

Release date 17 August 2026

Released by Nicole Geiger

Signature

pp. N. Geiger

REV 10_00 / 2023-09-15