

**Certificate of Analysis  
GM501 PVP**

Batchnumber	PVP-1815
Production date	20 August 2018
Expiry date	31 May 2019
Total volume (kits)	420
Batchnumber HSA	B3AKA00593

	<b>Results</b>	<b>Specifications</b>
pH	7,44	7,20-7,60
Osmolality	295	280-320 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,152	< 1,00 EU/mL
Viscosity in cP (26°C)	244,7	> 170 cP
Mouse Embryo Assay (1-cell)	95,24	≥ 80 %

*Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.*

Release date	04 September 2018
Released by	Dr. Julia Heinzmann
Signature	

*p.p. Julia Heinzmann*