

**Certificate of Analysis  
GM501 PVP**

Batchnumber	PVP-2006
Production date	23 March 2020
Expiry date	31 December 2020
Total volume (ml)	590
Batchnumber HSA	B3ALC00833

	<b>Results</b>	<b>Specifications</b>
pH	7,49	7,20-7,60
Osmolality	301	280-320 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 1,00 EU/mL
Viscosity in cP (26°C)	286,7	> 170 cP
Mouse Embryo Assay (1-cell)	80,95	≥ 80 %

*Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.*

Release date 07 April 2020

Released by Vanessa Bierlich

Signature