

**Certificate of Analysis  
GM501 PVP**

Batchnumber PVP-2303  
 Production date 15 February 2023  
 Expiry date 30 November 2023  
 Total volume (ml) 670  
 Batchnumber HSA B03G002753

	<b>Results</b>	<b>Specifications</b>
pH	7,51	7,20-7,60
Osmolality	303	280-320 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,5 EU/mL
Viscosity in cP (25°C)	251,9	> 170 cP
Mouse Embryo Assay (1-cell)	98	≥ 80 %

*Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.*

Release date 07 March 2023  
 Released by Nicole Geiger

Signature

