


Certificate of Analysis
GM501 PVP

Batchnumber	PVP-2311
Production date	08 June 2023
Expiry date	31 March 2024
Total volume (ml)	660
Batchnumber HSA	B03G034663

	Results	Specifications
pH	7,45	7,20-7,60
Osmolality	310	280-320 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,5 EU/mL
Viscosity in cP (25°C)	254,3	> 170 cP
Mouse Embryo Assay (1-cell)	90,48	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date	07 August 2023
Released by	Nicole Geiger
Signature	

REV 08_00 / 2022-12-01