

**Certificate of Analysis  
GM501 PVP**

Batchnumber	PVP-2316
Production date	08 August 2023
Expiry date	31 May 2024
Total volume (ml)	660
Batchnumber HSA	B03H003853

	<b>Results</b>	<b>Specifications</b>
pH	7,46	7,20-7,60
Osmolality	310	280-320 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,5 EU/mL
Viscosity in cP (25°C)	247,1	> 170 cP
Mouse Embryo Assay (1-cell)	85,71	≥ 80 %

*Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.*

Release date	23 August 2023
Released by	Nicole Geiger
Signature	

REV 07\_00 / 2022-12-01