

**Certificate of Analysis
GM501 PVP**

Batchnumber	PVP-2317
Production date	30 August 2023
Expiry date	31 May 2024
Total volume (ml)	660
Batchnumber HSA	B03H003853

	Results	Specifications
pH	7,45	7,20-7,60
Osmolality	302	280-320 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,5 EU/mL
Viscosity in cP (25°C)	198	> 170 cP
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date 14 September 2023

Released by Nicole Geiger

Signature

pp. N. Geiger

REV 07_00 / 2022-12-01