

**Certificate of Analysis
GM501 SpermMobil**

Batchnumber	SMOBIL-1904
Production date	11 June 2019
Expiry date	31 December 2019
Total volume (mL)	1780

	Results	Specifications
pH	7,37	7,20-7,60
Osmolality	298	280-320 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,50 EU/mL
Sperm survival Test 4 hours exposure	99	≥ 80 %

Release date 26 June 2019

Released by Vanessa Bierlich

Signature