

Certificate of Analysis
GM501 SpermMobil

Batchnumber	SMOBIL-2306
Production date	19 July 2023
Expiry date	31 January 2024
Total volume (mL)	2750

	Results	Specifications
pH	7,28	7,20-7,60
Osmolality	297	280-320 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,50 EU/mL
Sperm survival Test 4 hours exposure	100	≥ 80 %

Release date	04 August 2023
Released by	Nicole Geiger
Signature	

pp. N. Geiger