

**Certificate of Analysis
GM501 SpermMobil**

Batchnumber SMOBIL-2404
 Production date 03 June 2024
 Expiry date 31 December 2024
 Total volume (ml) 2680

Parameter	Results	Specifications
pH	7,31	7,20-7,60
Osmolality	298	280-320 mOsm/kg
Endotoxin level	< 0,05	< 0,50 EU/mL
Sterility test (according to the current Ph.Eur. 2.6.1 / USP <71>)	NO GROWTH	No growth
HSSA (% motility compared with control after 4h exposure to the test medium)	82	≥ 80 %

Release date 18 June 2024
 Released by Nicole Geiger
 Signature

pp. N. Geiger

REV 08_00 / 2023-09-15