

Certificate of Analysis GM501 SpermMobil

Batchnumber	SMOBIL-2502
Production date	04 March 2025
Expiry date	30 September 2025
Total volume (ml)	2300

Product was sterilized by aseptic processing techniques

Parameter	Results	Specifications
pH at 37°C	7,28	7,20-7,60
Osmolality	297	280-320 mOsm/kg
Endotoxin level	< 0,05	< 0,50 EU/mL
Sterility test (according to the current Ph.Eur. 2.6.1 / USP <71>)	No growth	No growth
HSSA (% motility compared with control after 4h exposure to the test medium)	90	≥ 80 %

Release date	19 March 2025
Released by	Joana Sturm
Signature	

REV 09_00 /2025-03-19