

Certificate of Analysis GM501 SpermMobil

Batchnumber	SMOBIL-2508
Production date	23 October 2025
Expiry date	30 April 2026
Total volume (ml)	2600
Product was sterilized by aseptic processing techniques	

Parameter	Results	Specifications
pH at 37°C	7,33	7,20-7,60
Osmolality	296	280-320 mOsm/kg
Endotoxin level	< 0,05	< 0,50 EU/mL
Sterility test (according to the current Ph.Eur. 2.6.1 / USP <71>)	No growth	No growth
HSSA (% motility compared with control after 4h exposure to the test medium)	99	≥ 80 %

Release date	07 November 2025
Released by	Nicole Geiger
Signature	

pp. P. Geiger

REV 09_00 /2025-03-19