

**Certificate of Analysis
GM501 SpermStore**

Batchnumber	SPSTORE-1913
Production date	11 June 2019
Expiry date	31 December 2020
Batchnumber HSA	K824B6871
Total volume (No. of bottles)	750

	Results	Specifications
pH	7,42	7,20-7,60
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Sperm survival Test 4 hours exposure	99	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date 26 June 2019

Released by Vanessa Bierlich

Signature

P.P. 