

**Certificate of Analysis  
GM501 SpermStore**

Batchnumber	SPSTORE-2114
Production date	09 August 2021
Expiry date	28 February 2023
Batchnumber HSA	B3ALE00033
Total volume (No. of bottles)	1160

	<b>Results</b>	<b>Specifications</b>
pH	7,52	7,20-7,60
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Sperm survival Test 4 hours exposure	99	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date 25 August 2021

Released by Dr. Julia Heinzmann

Signature

*p.p. Julia Heinzmann*