

**Certificate of Analysis  
GM501 SpermStore**

Batchnumber SPSTORE-2329  
 Production date 22 November 2023  
 Expiry date 31 May 2025  
 Batchnumber HSA B03H025323  
 Total number of kits 134

| Parameter                                                              | Result    | Specification |
|------------------------------------------------------------------------|-----------|---------------|
| pH                                                                     | 7,49      | 7,20-7,60     |
| Endotoxin level                                                        | < 0,1     | <0,25 EU/ml   |
| Concentration HSA                                                      | 3,9       | 3,2-4,8 g/l   |
| Sterility test<br>(according to the current Ph. Eur. 2.6.1 / USP <71>) | No growth | No growth     |
| Sperm Survival Test (4 hours exposure)                                 | 100       | ≥ 80 %        |

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date 07 December 2023

Released by Nicole Geiger

Signature



REV 07\_00 / 2023-09-15