


**Certificate of Analysis  
GM501 SpermStore**

Batchnumber SPSTORE-2418  
 Production date 19 July 2024  
 Expiry date 31 January 2026  
 Batchnumber HSA 297  
 Total number of kits B03H100413

Parameter	Result	Specification
pH	7,47	7,20-7,60
Endotoxin level	< 0,1	<0,25 EU/ml
Concentration HSA	3,9	3,2-4,8 g/l
Sterility test (according to the current Ph. Eur. 2.6.1 / USP <71>)	No growth	No growth
Sperm Survival Test (4 hours exposure)	96	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphillis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date 05 August 2024  
 Released by Nicole Geiger  
 Signature



REV 07\_00 / 2023-09-15