

Certificate of Analysis GM501 VitriStore Thaw

Batchnumber	VITRI-1910
Production date	12 August 2019
Expiry date	31 August 2020
Batchnumber HSA	B3ALC00833
Total volume (Kits)	199

	Results	Specifications
GM501 VitriStore Thaw 1		
pH	7,32	7,20-7,40
Sterility test	PASSED	PASSED
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	835	805-850 mOsm/kg

	Results	Specifications
GM501 VitriStore Thaw 2		
pH	7,34	7,20-7,40
Sterility test	PASSED	PASSED
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	554	535-565 mOsm/kg

	Results	Specifications
GM501 VitriStore Thaw 3		
pH	7,28	7,20-7,40
Sterility test	PASSED	PASSED
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	413	405-435 mOsm/kg

	Results	Specifications
GM501 VitriStore Thaw 4		
pH	7,32	7,20-7,40
Sterility test	PASSED	PASSED
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	282	270-290 mOsm/kg

Mouse Embryo Assay (1-cell)	100	≥ 80 %
-----------------------------	-----	--------

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date	30 August 2019
--------------	----------------

Released by	Dr. Julia Heinzmann
-------------	---------------------

Signature

