

## Certificate of Analysis GM501 VitriStore Thaw

Batchnumber	VITRI-2006
Production date	26 May 2020
Expiry date	31 May 2021
Batchnumber HSA	B3ALD01343
Total volume (Kits)	199

	Results	Specifications
<b>GM501 VitriStore Thaw 1</b>		
pH	7,36	7,20-7,40
Sterility test	PASSED	PASSED
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	822	805-850 mOsm/kg

	Results	Specifications
<b>GM501 VitriStore Thaw 2</b>		
pH	7,35	7,20-7,40
Sterility test	PASSED	PASSED
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	544	535-565 mOsm/kg

	Results	Specifications
<b>GM501 VitriStore Thaw 3</b>		
pH	7,35	7,20-7,40
Sterility test	PASSED	PASSED
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	418	405-435 mOsm/kg

	Results	Specifications
<b>GM501 VitriStore Thaw 4</b>		
pH	7,36	7,20-7,40
Sterility test	PASSED	PASSED
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	282	270-290 mOsm/kg

Mouse Embryo Assay (1-cell)	100	≥ 80 %
-----------------------------	-----	--------

*Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.*

Release date	16 June 2020
--------------	--------------

Released by	Dr. Julia Heinzmann
-------------	---------------------

Signature

*p.p. Julia Heinzmann*