

Certificate of Analysis GM501 VitriStore Thaw

Batchnumber	VITRI-2008
Production date	06 August 2020
Expiry date	31 August 2021
Batchnumber HSA	B3ALD01343
Total volume (Kits)	200

	Results	Specifications
GM501 VitriStore Thaw 1		
pH	7,32	7,20-7,40
Sterility test	PASSED	PASSED
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	834	805-850 mOsm/kg

	Results	Specifications
GM501 VitriStore Thaw 2		
pH	7,31	7,20-7,40
Sterility test	PASSED	PASSED
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	552	535-565 mOsm/kg

	Results	Specifications
GM501 VitriStore Thaw 3		
pH	7,31	7,20-7,40
Sterility test	PASSED	PASSED
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	413	405-435 mOsm/kg

	Results	Specifications
GM501 VitriStore Thaw 4		
pH	7,31	7,20-7,40
Sterility test	PASSED	PASSED
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	278	270-290 mOsm/kg

Mouse Embryo Assay (1-cell)	100	≥ 80 %
-----------------------------	-----	--------

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date	28 August 2020
--------------	----------------

Released by	Vanessa Bierlich
-------------	------------------

Signature

