

Certificate of Analysis

GM501 VitriStore Thaw

Batchnumber	VITRI-2319
Production date	23 August 2023
Expiry date	31 August 2024
Batchnumber HSA	B03H003853
Total volume (Kits)	172

GM501 VitriStore Thaw 1	Results	Specifications
pH	7,32	7,20-7,40
Sterility test	PASSED	PASSED
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	838	805-850 mOsm/kg

GM501 VitriStore Thaw 2	Results	Specifications
pH	7,28	7,20-7,40
Sterility test	PASSED	PASSED
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	557	535-565 mOsm/kg

GM501 VitriStore Thaw 3	Results	Specifications
pH	7,29	7,20-7,40
Sterility test	PASSED	PASSED
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	419	405-435 mOsm/kg

GM501 VitriStore Thaw 4	Results	Specifications
pH	7,3	7,20-7,40
Sterility test	PASSED	PASSED
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	285	270-290 mOsm/kg

Mouse Embryo Assay (1-cell)	100	≥ 80 %
-----------------------------	-----	--------

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date	12 September 2023
--------------	-------------------

Released by	Nicole Geiger
-------------	---------------

Signature

pp. N. Geiger