

Certificate of Analysis GM501 VitriStore Thaw

Batchnumber	VITRI-2613
Production date	10 June 2026
Expiry date	30 June 2027
Total batch size	210 Kits
Batchnumber HSA	B03J125033

Parameter	Results	Specifications
GM501 VitriStore Thaw 1		
pH	7,33	7,20-7,40
Sterility test ¹	No growth	No growth
Endotoxin	< 0,15	< 0,25 EU/mL
Osmolality	827	805-850 mOsm/kg
Concentration HSA	17,8	14,5-21,7 g/l
GM501 VitriStore Thaw 2		
pH	7,31	7,20-7,40
Sterility test ¹	No growth	No growth
Endotoxin	< 0,15	< 0,25 EU/mL
Osmolality	555	535-565 mOsm/kg
Concentration HSA	18,7	15,2-22,8 g/l
GM501 VitriStore Thaw 3		
pH	7,26	7,20-7,40
Sterility test ¹	No growth	No growth
Endotoxin	< 0,15	< 0,25 EU/mL
Osmolality	421	405-435 mOsm/kg
Concentration HSA	20,1	15,5-23,3 g/l
GM501 VitriStore Thaw 4		
pH	7,31	7,20-7,40
Sterility test ¹	No growth	No growth
Endotoxin	< 0,15	< 0,25 EU/mL
Osmolality	285	270-290 mOsm/kg
Concentration HSA	20,5	16,0-24,0 g/l

Mouse Embryo Assay (1-cell)	86	≥ 80 %
-----------------------------	----	--------


Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

¹(according to the current Ph.Eur.2.6.1 / USP <71>)

Release date	26 June 2026
--------------	--------------

Released by	Dr. René Richarz
-------------	------------------

Signature



REV 06_00 / 2023-09-15