

Certificate of Analysis

GM501 WASH with Phenol Red and Gentamicin

Batchnumber	WASH-1808PRG
Production date	26 July 2018
Expiry date	31 January 2019
Total volume (mL)	35500
Batchnumber HSA	M749A6671 / K809B6671

	Results	Specifications
pH	7,24	7,20-7,50
Osmolality	280	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	90,48	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date	10.08.0218
--------------	------------

Released by	Dr. Julia Heinzmann
-------------	---------------------

Signature

p.p. Julia Heinzmann