

**Certificate of Analysis  
GM501 WASH with Phenol Red and Gentamicin**

Batchnumber	WASH-2109PRG
Production date	07 July 2021
Expiry date	31 January 2022
Batchnumber HSA	33500
Total volume (mL)	B3ALE00033

	<b>Results</b>	<b>Specifications</b>
pH	7,22	7,20-7,50
Osmolality	280	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

*Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.*

Release date 22 July 2021

Released by Dr. Julia Heinzmann

Signature

*p.p. Julia Heinzmann*