

**Certificate of Analysis**  
**GM501 WASH**

Batchnumber	WASH-2307
Production date	23 May 2023
Expiry date	30 November 2023
Batchnumber HSA	B03G034663
Total volume (mL)	19800

	<b>Results</b>	<b>Specifications</b>
pH	7,31	7,20-7,50
Osmolality	280	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

*Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.*

Release date	13 June 2023
Released by	Nicole Geiger
Signature	

pp. N. Geiger