

**Certificate of Analysis**  
**GM501 SpermActive**

Batchnumber	WASH-2308PRG
Production date	05 July 2023
Expiry date	31 January 2024
Batchnumber HSA	B03H003853
Total volume (mL)	32200

	<b>Results</b>	<b>Specifications</b>
pH	7,34	7,20-7,50
Osmolality	282	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	95	≥ 80 %

*Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.*

Release date	25 July 2023
Released by	Dr. Julia Heinzmann
Signature	

*p.p. Julia Heinzmann*

REV 05\_00 / 2019-01-15