

## Certificate of Analysis GM501 WASH

Batchnumber	WASH-2403
Production date	06 May 2024
Expiry date	30 November 2024
Batchnumber HSA	B03H100413
Total volume (mL)	10200

Parameter	Results	Specifications
pH	7,34	7,20-7,50
Osmolality	280	270-290 mOsm/kg
Bacterial screen	No growth	No growth
Fungal screen	No growth	No growth
Endotoxin	< 0,05	< 0,25 EU/mL
HSA concentration	4,9	4,0-6,0 mg/ml
Mouse Embryo Assay (1-cell)	100	≥ 80 %

*Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.*

Release date	21 May 2024
Released by	Joana Sturm
Signature	

REV 07\_00 / 2024-02-08