

**Certificate of Analysis
GM501 WASH with Phenol Red and Gentamicin**

Batchnumber	WASH-2403PRG
Production date	06 May 2024
Expiry date	30 November 2024
Batchnumber HSA	B03H100413
Total volume (mL)	33400

Parameter	Results	Specifications
pH	7,32	7,20-7,50
Osmolality	281	270-290 mOsm/kg
Bacterial screen	No growth	No growth
Fungal screen	No growth	No growth
Endotoxin	< 0,05	< 0,25 EU/mL
HSA concentration	4,9	4,0-6,0 mg/ml
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date 21 May 2024

Released by Joana Sturm

Signature



REV 07_00 / 2024-02-08