

**Certificate of Analysis
GM501 SpermActive**

Batchnumber WASH-2405PRG
 Production date 19 July 2024
 Expiry date 31 January 2025
 Batchnumber HSA B03H100413
 Total volume (mL)


Parameter	Results	Specifications
pH	7,31	7,20-7,50
Osmolality	278	270-290 mOsm/kg
Bacterial screen	No growth	No growth
Fungal screen	No growth	No growth
Endotoxin	< 0,05	< 0,25 EU/mL
HSA concentration	4,7	4,0-6,0 mg/ml
Mouse Embryo Assay (1-cell)	95,24	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date 05 August 2024

Released by Nicole Geiger

Signature



REV 07_00 / 2024-02-08