

**Certificate of Analysis**  
**GM501 WASH with Phenol Red and Gentamicin**


Batchnumber	WASH-2406PRG
Production date	02 September 2024
Expiry date	31 March 2025
Batchnumber HSA	B03H100413
Total volume (mL)	27600

<b>Parameter</b>	<b>Results</b>	<b>Specifications</b>
pH	7,26	7,20-7,50
Osmolality	281	270-290 mOsm/kg
Endotoxin level	< 0,05	< 0,25 EU/ml
Concentration HSA	4,7	4,0-6,0 g/l
Concentration Gentamicin	12	5-15 mg/l
Sterility test (performed according to the current Ph.Eur.2.6.1 / USP <71>)	No growth	No growth
Mouse Embryo Assay (1-cell to expanded blastocyst within 96 hours)	95	≥ 80 %

*Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.*

Release date	18 September 2024
Released by	Nicole Geiger

Signature



REV 08\_00 / 2024-09-17