

**Certificate of Analysis  
GM501 WASH**

Batchnumber WASH-2505  
 Production date 20 June 2025  
 Expiry date 31 December 2025  
 Total volume (ml) 5300  
 Batchnumber HSA B03H121903

Product was sterilized by aseptic processing techniques

Parameter	Results	Specifications
pH (@37°C - 6% CO <sub>2</sub> )	7,31	7,20-7,50
Osmolality	280	270-290 mOsm/kg
Endotoxin level	< 0,05	< 0,25 EU/ml
Concentration HSA	4,8	4,0-6,0 g/l
Sterility test (performed according to the current Ph.Eur.2.6.1 / USP <71>)	No growth	No growth
Mouse Embryo Assay (1-cell to expanded blastocyst within 96 hours)	100	≥ 80 %

*Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.*

Release date 08 July 2025  
 Released by Nicole Geiger  
 Signature

pp. N. Geiger

REV 09\_00 / 2025-03-19