

Certificate of Analysis GM501 SpermActive


Batchnumber	WASH-2602PRG
Production date	12 February 2026
Expiry date	31 August 2026
Total volume (mL)	33700
Batchnumber HSA	B03J103333

Product was sterilized by aseptic processing techniques ☑

Parameter	Results	Specifications
pH (@37°C - 6% CO ₂)	7,28	7,20-7,50
Osmolality	279	270-290 mOsm/kg
Endotoxin level	< 0,05	< 0,25 EU/ml
Concentration HSA	4,6	4,0-6,0 g/l
Concentration gentamicin	10	5-15 mg/l
Sterility test (performed according to the current Ph.Eur.2.6.1 / USP <71>)	No growth	No growth
Mouse Embryo Assay (1-cell to expanded blastocyst within 96 hours)	100	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date	27 February 2026
Released by	Nicole Geiger
Signature	



REV 09_00 / 2025-03-19