

**Certificate of Analysis**  
**GM501 WASH with Phenol Red and Gentamicin**

Batchnumber WASH-2603PRG  
 Production date 10 March 2026  
 Expiry date 30 September 2026  
 Total volume (mL) 33100  
 Batchnumber HSA B03J103333

Product was sterilized by aseptic processing techniques

<b>Parameter</b>	<b>Results</b>	<b>Specifications</b>
pH (@37°C - 6% CO <sub>2</sub> )	7,26	7,20-7,50
Osmolality	280	270-290 mOsm/kg
Endotoxin level	< 0,05	< 0,25 EU/ml
Concentration HSA	4,7	4,0-6,0 g/l
Concentration gentamicin	9	5-15 mg/l
Sterility test (performed according to the current Ph.Eur.2.6.1 / USP <71>)	No growth	No growth
Mouse Embryo Assay (1-cell to expanded blastocyst within 96 hours)	100	≥ 80 %

*Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.*

Release date 25 March 2026  
 Released by Dr. René Richarz  
 Signature



REV 09\_00 / 2025-03-19