

Certificate of Analysis

GM501 WASH with Phenol Red and Gentamicin

Batchnumber WASH-2609PRG

Production date 16 September 2026

Expiry date 31 March 2027

Total volume (mL) 22600

Batchnumber HSA B03J125033

Product was sterilized by aseptic processing techniques


Parameter	Results	Specifications
pH (@37°C - 6% CO ₂)	7,25	7,20-7,50
Osmolality	280	270-290 mOsm/kg
Endotoxin level	< 0,05	< 0,25 EU/ml
Concentration HSA	4,7	4,0-6,0 g/l
Concentration gentamicin	8	5-15 mg/l
Sterility test (performed according to the current Ph.Eur.2.6.1 / USP <71>)	No growth	No growth
Mouse Embryo Assay (1-cell to expanded blastocyst within 96 hours)	100	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date 02 October 2026

Released by Dr. René Richarz

Signature



REV 09_00 / 2025-03-19