

**Certificate of Analysis
GM501 SpermMobil**

Batchnumber	SMOBIL-1707
Production date	01 December 2017
Expiry date	30 June 2018
Total volume (mL)	1710

	Results	Specifications
pH	7,33	7,20-7,60
Osmolality	300	280-320 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,50 EU/mL
Sperm survival Test 4 hours exposure	88	≥ 80 %

Release date 18 December 2017

Released by Dr. Julia Heinzmann

Signature

p.p. Julia Heinzmann