



SparMED Aps  
 Toppevadvej 34-38  
 DK-3660  
 Stenlose, Denmark



140 Hale Street  
 Haverhill, MA 01830  
[qc@embryotech.com](mailto:qc@embryotech.com)

ELI Accession Number: SPAR-4065-1215

Date of completion: 12-22-2015

Lot numbers: 07603

Reference numbers: OOPW-SC01, OOPW-CT01  
 OOPW-AT10, OOPW-OT10

**Description of test article(s):**

Oosafe® Sperm Collection Cup, Centrifuge Tube, Andrology Tube, OPU Tube

**Assay system requested by customer:** 1mL of culture medium was placed in each of the test articles (3) (OOPW-SC01, OOPW-CT01, OOPW-AT10) and incubated at 37°C for 30-minutes. Post incubation the culture medium was extracted from each test article and pooled. 1mL of the extracted culture medium was expelled into the fourth test article OOPW-OT10 overlaid with oil. One-cell mouse embryos were then added to each well of the test article and cultured for 96-hours.

**Control assay method and results:** 15 1-cell (B6C3F1 X B6D2F1) embryos were cultured in 0.5mL drops in a 60x15mm dish using culture medium:

15 / 15 (100 %)  
 15 / 15 (100 %)

1-cell to 2-cell within 24 hr  
 1-cell to expanded blastocyst within 96 hr

*For a valid assay, Embryotech™ requires at least 70% of 1-cell stage control embryos to develop to expanded blastocyst within 96-hours.*

**Test assay method and results:** 21 1-cell (B6C3F1 X B6D2F1) embryos were cultured in one of the test articles using the extracted culture medium:

21 / 21 (100 %)  
 21 / 21 (100 %)

1-cell to 2-cell within 24 hr  
 1-cell to expanded blastocyst within 96 hr

**Summary of observations:** All test and control embryos were selected randomly from a common pool of freshly collected embryos and were cultured in the same incubator at 37°C and 5.0% CO<sub>2</sub>. 100 percent of the control embryos developed to the expanded blastocyst stage within 96-hours. 10 percent of the embryos cultured in the test article developed to the expanded blastocyst stage within 96-hours.

signature  
 Study Director

12-24-2015  
 date

signature  
 Quality Reviewer

12-24-2015  
 date