

## Certificate of Analysis

**Date of issue:** 23 Jan 2019

**Product ID:** Oosafe® Plasticware: OOPW-SC01

**LOT No.:** 07981

**Expiry date:** 2023-12

**Storage Conditions:** 20°C – 30°C, dry room, no exposal to sun-light

**Quality Assurance:**

| Analyses   | Results |
|--|---------|
| Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours. | Passed  |
| Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.   | Passed  |
| Proven non-toxic by Limulus Amebocyte Lysate (LAL) test. Pass criteria <0.03 EU/device.  | Passed  |
| Sterilization by gamma irradiation. Delivered irradiation dose: 8.6 kGy-9.5 kGy. Specified irradiation dose: 8.0 kGy-10.0 kGy.   | Passed  |

Quality control according to the ISO 13485:2012

**GOosafe with SparMED!**

**Date:** 23 January 2019

Simona Laurinavičiute  
Quality Control Specialist  
SparMED ApS

