

Certificate of Analysis

Date of issue: 23 Jan 2019

Product ID: Oosafe® Plasticware: OOPW-SW03

LOT No.: 07982

Expiry date: 2023-12

Storage Conditions: 20°C – 30°C, dry room, no exposal to sun-light

Quality Assurance:

Analyses	Results
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	Passed
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.	Passed
Proven non-toxic by Limulus Amebocyte Lysate (LAL) test. Pass criteria <0.03 EU/device.	Passed
Sterilization by gamma irradiation. Delivered irradiation dose: 8.6 kGy-9.5 kGy. Specified irradiation dose: 8.0 kGy-10.0 kGy.	Passed

Quality control according to the ISO 13485:2012

GOosafe with SparMED!

Date: 23 January 2019

Simona Laurinavičiute
Quality Control Specialist
SparMED ApS

