

Certificate of Analysis

Product ID: OOSAFE® 35 MM DISH, HIGH WALL

Product code: OOPW-TF03

LOT No.: 07993

Expiry date: 2024-04

Storage Conditions: 20-30°C, dry room, no exposure to sunlight

Quality Assurance:

| Analysis | Reference No. | Results | Completion date |
|----------------------------------------------------------------------------------------------------------------------------------|----------------|---------|-----------------|
| Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours. | SPAR-2192-0719 | Passed | 26-07-2019 |
| Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven. | S3964-0719SPAR | Passed | 24-07-2019 |
| Proven non-toxic by Limulus Amebocyte Lysate (LAL) test. Pass criteria <0.03 EU/device. | E9467-0719SPAR | Passed | 22-07-2019 |
| Sterilization by gamma irradiation. Delivered irradiation dose: 8.6 kGy-9.5 kGy. Specified irradiation dose: 8.0 kGy-10.0 kGy. | N/A | Passed | N/A |

Quality Control according to ISO 13485:2016

GOosafe with SparMED!

Inga Sileikaite,
Quality Control Department
SparMED ApS



Date: 29-07-2019

