

Certificate of Analysis

Product ID: OOSAFE® 4 WELL DISH NON-TREATED SURFACE

Product code: OOPW-FW04

LOT No.: 08015

Expiry date: 2024-09

Storage Conditions: 20-30°C, dry room, no exposure to sunlight

Quality Assurance:

Analysis	Reference No.	Results	Completion date
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR-2642-1019	Passed	07-10-2019
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.	S4071-1019SPAR	Passed	05-10-2019
Proven non-toxic by Limulus Amebocyte Lysate (LAL) test. Pass criteria <0.03 EU/device.	E9646-1019SPAR	Passed	03-10-2019
Sterilization by gamma irradiation. Delivered irradiation dose: 8.6 kGy-9.5 kGy. Specified irradiation dose: 8.0 kGy-10.0 kGy.	NA	Passed	NA

Quality Control according to ISO 13485:2016

GOosafe with SparMED!

Inga Sileikaite,
Quality Control Department
SparMED ApS

Date: 08-10-2019

