

## Certificate of Analysis

**Product ID:** Oosafe® ICSI/IMSI Dish for Sperm Selection

**Product code:** OOPW-IC03

**LOT No.:** 08055

**Expiry date:** 2025-04

**Storage Conditions:** 20-30°C, dry room, no exposure to sunlight

**Quality Assurance:**

Analysis	Reference No.	Results	Completion date
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR-4079-0620	Passed	06-06-2020
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.	S4478-0620SPAR	Passed	05-06-2020
Proven non-toxic by Limulus Amebocyte Lysate (LAL) test. Pass criteria <0.03 EU/device.	E0252-0620SPAR	Passed	02-06-2020
Sterilization by gamma irradiation. Delivered irradiation dose: 8.6 kGy-9.5 kGy. Specified irradiation dose: 8.0 kGy-10.0 kGy.	N/A	Passed	N/A

Quality Control according to ISO 13485:2016

### GOosafe with SparMED!

Florencia Garcia  
Quality Control Department  
SparMED ApS  
**Date:** 03-08-2020

