

## Certificate of Analysis

**Product ID:** Oosafe® OPU Tube, 14 ml**Product code:** OOPW-OT10**LOT No.:** 08095**Expiry date:** 2025-07**Storage Conditions:** 20-30°C, dry room, no exposure to sunlight**Quality Assurance:**

Analysis	Reference No.	Results	Completion date
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR-4461-0920	Passed	08-09-2020
Proved stable human sperm motility: $\geq 75\%$ sperm motility after 24 hours proven.	S4596-0920SPAR	Passed	09-09-2020
Proven non-toxic by Limulus Amebocyte Lysate (LAL) test. Pass criteria $< 0.03$ EU/device.	E0441-0920SPAR	Passed	09-09-2020
Sterilization by gamma irradiation. Delivered irradiation dose: 8.6 kGy-9.5 kGy. Specified irradiation dose: 8.0 kGy-10.0 kGy.	N/A	Passed	N/A

Quality Control according to ISO 13485:2016

**GOosafe with SparMED!**

Florencia Garcia  
Quality Control Department  
SparMED ApS  
**Date:** 10-09-2020

