

## Certificate of Analysis

**Product ID: OOSAFE® OPU Tube, 14 ml**

**Product code: OOPW-OT10**

**LOT No.: 08213**

**Expiry date: 2027-01**

**Storage Conditions: 20-30°C, dry room, no exposure to sunlight**

**Quality Assurance:**

Analysis	Reference No.	Results	Completion date
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR-7568-0122	Passed	04-02-2022
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.	S5621-0122SPAR	Passed	10-02-2022
Proven non-toxic by Limulus Amebocyte Lysate (LAL) test. Pass criteria <0.03 EU/device.	E1940-0122SPAR	Passed	01-02-2022
Sterilization by gamma irradiation. Delivered irradiation dose: 8.6 kGy-9.5 kGy. Specified irradiation dose: 8.0 kGy-10.0 kGy.	N/A	Passed	N/A

Quality Control according to ISO 13485:2016

### GOosafe with SparMED!

Zahraa Baker  
Quality Control Department  
SparMED ApS  
Date: 14-02-2022

